

DAB2 Polyclonal Antibody FITC Conjugated

Catalog No: #C92416FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DAB2 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	31-130/770
Conjugates	FITC
Target Name	DAB2
Other Names	DAB2
Accession No.	Uniprot:P98082Gene ID:1601
Uniprot	P98082
GeneID	1601
Calculated MW	82
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only