

## DAL1/EPB41L3 Polyclonal Antibody FITC Conjugated

Catalog No: #C92417FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DAL1/EPB41L3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	131-230/1087
Conjugates	FITC
Target Name	DAL1/EPB41L3
Other Names	DAL1/EPB41L3
Accession No.	Uniprot:Q9Y2J2Gene ID:23136
Uniprot	Q9Y2J2
GeneID	23136
Calculated MW	120
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only