

TSLC1 Polyclonal Antibody PE Conjugated

Catalog No: #C92418PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TSLC1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	65-160/442
Conjugates	PE
Target Name	TSLC1
Other Names	TSLC1
Accession No.	Uniprot:Q9BY67Gene ID:23705
Uniprot	Q9BY67
GeneID	23705
Calculated MW	44
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only