

## CD16 Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C92420PE-Cy7

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD16 Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	131-230/254
Conjugates	PE-Cy7
Target Name	CD16
Other Names	CD16
Accession No.	Uniprot:P08637Gene ID:2214
Uniprot	P08637
GeneID	2214
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only