

## RBMX/hnRNP G Polyclonal Antibody PE Conjugated

Catalog No: #C92431PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RBMX/hnRNP G Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/391
Conjugates	PE
Target Name	RBMX/hnRNP G
Other Names	RBMX/hnRNP G
Accession No.	Uniprot:P38159Gene ID:27316
Uniprot	P38159
GeneID	27316
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only