

TM4SF3 Polyclonal Antibody FITC Conjugated

Catalog No: #C92450FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TM4SF3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/237
Conjugates	FITC
Target Name	TM4SF3
Other Names	TM4SF3
Accession No.	Uniprot:P19075Gene ID:7103
Uniprot	P19075
GeneID	7103
Calculated MW	26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only