

Prostaglandin E Receptor EP1 Polyclonal Antibody PE Conjugated



Catalog No: #C92459PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostaglandin E Receptor EP1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	61-160/402
Conjugates	PE
Target Name	Prostaglandin E Receptor EP1
Other Names	Prostaglandin E Receptor EP1
Accession No.	Uniprot:P34995Gene ID:5731
Uniprot	P34995
GeneID	5731
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only