

## ICAM1 Polyclonal Antibody HRP Conjugated

Catalog No: #C92462HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ICAM1 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu Dog Pig
Immunogen Type	351-450/532
Conjugates	HRP
Target Name	ICAM1
Other Names	ICAM1
Accession No.	Uniprot:P05362Gene ID:3383
Uniprot	P05362
GeneID	3383
Calculated MW	55
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only