

Filaggrin Polyclonal Antibody FITC Conjugated

Catalog No: #C92463FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Filaggrin Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-150/4061
Conjugates	FITC
Target Name	Filaggrin
Other Names	Filaggrin
Accession No.	Uniprot:P20930Gene ID:2312
Uniprot	P20930
GeneID	2312
Calculated MW	447
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only