

SCGB3A1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92468Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SCGB3A1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu
Immunogen Type	51-104/104
Conjugates	Bio
Target Name	SCGB3A1
Other Names	SCGB3A1
Accession No.	Uniprot:Q96QR1Gene ID:92304
Uniprot	Q96QR1
GeneID	92304
Calculated MW	8.2
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only