

SMARCA2 Polyclonal Antibody APC Conjugated

Catalog No: #C92469APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | SMARCA2 Polyclonal Antibody APC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IF |
| Species Reactivity | Hu Ms Ra |
| Immunogen Type | 1401-1590/1590 |
| Conjugates | APC |
| Target Name | SMARCA2 |
| Other Names | SMARCA2 |
| Accession No. | Uniprot:P51531Gene ID:6595 |
| Uniprot | P51531 |
| GeneID | 6595 |
| Calculated MW | 230 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Note: This product is for in vitro research use only