

## HIV2 gp36 Polyclonal Antibody HRP Conjugated

Catalog No: #C92470HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	HIV2 gp36 Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB,IHC
Species Reactivity	HIV2
Immunogen Type	78-125/858
Conjugates	HRP
Target Name	HIV2 gp36
Other Names	HIV2 gp36
Accession No.	Uniprot:P04577Gene ID:
Uniprot	P04577
Calculated MW	43/92
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:500-2000 IHC-F=1:50-200

Note: This product is for in vitro research use only