

CDH11/OB Cadherin Polyclonal Antibody Biotin Conjugated

Catalog No: #C92477Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CDH11/OB Cadherin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-300/796
Conjugates	Bio
Target Name	CDH11/OB Cadherin
Other Names	CDH11/OB Cadherin
Accession No.	Uniprot:P55287Gene ID:1009
Uniprot	P55287
GeneID	1009
Calculated MW	90
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only