

## PRLR Polyclonal Antibody APC Conjugated

Catalog No: #C92478APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PRLR Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	251-350/622
Conjugates	APC
Target Name	PRLR
Other Names	PRLR
Accession No.	Uniprot:P16471Gene ID:5618
Uniprot	P16471
GeneID	5618
Calculated MW	67
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only