

ITPR3 Polyclonal Antibody PE Conjugated

Catalog No: #C92487PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ITPR3 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/2671
Conjugates	PE
Target Name	ITPR3
Other Names	ITPR3
Accession No.	Uniprot:Q14643Gene ID:3708
Uniprot	Q14643
GeneID	3708
Calculated MW	294
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only