

# GABA A Receptor alpha 3 Polyclonal Antibody FITC Conjugated



Catalog No: #C92489FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GABA A Receptor alpha 3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/492
Conjugates	FITC
Target Name	GABA A Receptor alpha 3
Other Names	GABA A Receptor alpha 3
Accession No.	Uniprot:P34903Gene ID:2556
Uniprot	P34903
GeneID	2556
Calculated MW	52
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only