

Sema3C Polyclonal Antibody FITC Conjugated

Catalog No: #C92517FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sema3C Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/751
Conjugates	FITC
Target Name	Sema3C
Other Names	Sema3C
Accession No.	Uniprot:Q99985Gene ID:10512
Uniprot	Q99985
GeneID	10512
Calculated MW	80
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only