

SOX11 Polyclonal Antibody FITC Conjugated

Catalog No: #C92520FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SOX11 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/441
Conjugates	FITC
Target Name	SOX11
Other Names	SOX11
Accession No.	Uniprot:P35716Gene ID:6664
Uniprot	P35716
GeneID	6664
Calculated MW	47
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only