

NPC1 Polyclonal Antibody FITC Conjugated

Catalog No: #C92525FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | NPC1 Polyclonal Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 1181-1278/1287 |
| Conjugates | FITC |
| Target Name | NPC1 |
| Other Names | NPC1 |
| Accession No. | Uniprot:O15118Gene ID:4864 |
| Uniprot | O15118 |
| GeneID | 4864 |
| Calculated MW | 138 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only