

CCBR1 Polyclonal Antibody PE Conjugated

Catalog No: #C92542PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CCBR1 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/501
Conjugates	PE
Target Name	CCBR1
Other Names	CCBR1
Accession No.	Uniprot:Q9UPY5Gene ID:23657
Uniprot	Q9UPY5
GeneID	23657
Calculated MW	55
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only