

## BMP9 Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C92545PE-Cy5.5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	BMP9 Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	351-429/429
Conjugates	PE-Cy5.5
Target Name	BMP9
Other Names	BMP9
Accession No.	Uniprot:Q9UK05Gene ID:2658
Uniprot	Q9UK05
GeneID	2658
Calculated MW	12/47
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only