

phospho-PKC Epsilon(Ser729) Polyclonal Antibody HRP Conjugated

Catalog No: #C92553HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	phospho-PKC Epsilon(Ser729) Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC
Species Reactivity	Hu Ms Rat
Immunogen Type	GF(p-S)YF
Conjugates	HRP
Target Name	phospho-PKC Epsilon(Ser729)
Other Names	phospho-PKC Epsilon α 1/2Ser729 α 1/2
Accession No.	Uniprot:Q02156Gene ID:5581
Uniprot	Q02156
GeneID	5581
Calculated MW	84
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only