

RNA polymerase II Polyclonal Antibody PE Conjugated



Catalog No: #C92557PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	RNA polymerase II Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/1970
Conjugates	PE
Target Name	RNA polymerase II
Other Names	RNA polymerase II
Accession No.	Uniprot:P24928Gene ID:5430
Uniprot	P24928
GeneID	5430
Calculated MW	217
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only