

Neutrophil Elastase Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C92562PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neutrophil Elastase Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/267
Conjugates	PE-Cy7
Target Name	Neutrophil Elastase
Other Names	Neutrophil Elastase
Accession No.	Uniprot:P08246Gene ID:1991
Uniprot	P08246
GeneID	1991
Calculated MW	26
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only