

Wilms Tumor Protein Polyclonal Antibody FITC Conjugated

Catalog No: #C92563FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Wilms Tumor Protein Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-400/449
Conjugates	FITC
Target Name	Wilms Tumor Protein
Other Names	Wilms Tumor Protein
Accession No.	Uniprot:P19544Gene ID:7490
Uniprot	P19544
GeneID	7490
Calculated MW	55
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only