

EphA3 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92570Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EphA3 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	401-500/983
Conjugates	Bio
Target Name	EphA3
Other Names	EphA3
Accession No.	Uniprot:P29320Gene ID:2042
Uniprot	P29320
GeneID	2042
Calculated MW	106
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only