

EphA7 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92571Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EphA7 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	181-280/998
Conjugates	Bio
Target Name	EphA7
Other Names	EphA7
Accession No.	Uniprot:Q15375Gene ID:2045
Uniprot	Q15375
GeneID	2045
Calculated MW	50/107
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only