

EMR1 Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C92576PE-Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EMR1 Polyclonal Antibody PE-Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	701-800/886
Conjugates	PE-Cy5
Target Name	EMR1
Other Names	EMR1
Accession No.	Uniprot:Q14246Gene ID:2015
Uniprot	Q14246
GeneID	2015
Calculated MW	95
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only