

## MUC5AC Polyclonal Antibody Biotin Conjugated

Catalog No: #C92585Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MUC5AC Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	601-700/5030
Conjugates	Bio
Target Name	MUC5AC
Other Names	MUC5AC
Accession No.	Uniprot:P98088Gene ID:4586
Uniprot	P98088
GeneID	4586
Calculated MW	523
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only