

## CD200R Polyclonal Antibody APC Conjugated

Catalog No: #C92591APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD200R Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu
Immunogen Type	151-250/325
Conjugates	APC
Target Name	CD200R
Other Names	CD200R
Accession No.	Uniprot:Q6Q8B3Gene ID:131450
Uniprot	Q6Q8B3
GeneID	131450
Calculated MW	34
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only