

# Topoisomerase III Polyclonal Antibody FITC Conjugated

Catalog No: #C92592FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Topoisomerase III Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	271-370/1001
Conjugates	FITC
Target Name	Topoisomerase III
Other Names	Topoisomerase III
Accession No.	Uniprot:Q13472Gene ID:7156
Uniprot	Q13472
GeneID	7156
Calculated MW	112
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only