

beta Defensin 3 Polyclonal Antibody FITC Conjugated

Catalog No: #C92594FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	beta Defensin 3 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu
Immunogen Type	23-67/67
Conjugates	FITC
Target Name	beta Defensin 3
Other Names	beta Defensin 3
Accession No.	Uniprot:P81534Gene ID:414325
Uniprot	P81534
GeneID	414325
Calculated MW	5
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only