

## stabilin1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92601Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	stabilin1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	1201-1400/2570
Conjugates	Bio
Target Name	stabilin1
Other Names	stabilin1
Accession No.	Uniprot:Q9NY15Gene ID:23166
Uniprot	Q9NY15
GeneID	23166
Calculated MW	273
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200

Note: This product is for in vitro research use only