

Insulin Receptor R Polyclonal Antibody Biotin Conjugated

Catalog No: #C92606Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin Receptor R Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	801-900/1297
Conjugates	Bio
Target Name	Insulin Receptor R
Other Names	Insulin Receptor R
Accession No.	Uniprot:P14616Gene ID:3645
Uniprot	P14616
GeneID	3645
Calculated MW	129
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only