

# PGC1 alpha / beta Polyclonal Antibody APC Conjugated

Catalog No: #C92608APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PGC1 alpha / beta Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	151-250/798
Conjugates	APC
Target Name	PGC1 alpha / beta
Other Names	PGC1 alpha / beta
Accession No.	Uniprot:Q9UBK2Gene ID:10891
Uniprot	Q9UBK2
GeneID	10891
Calculated MW	88/113
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only