

VDAC3 Polyclonal Antibody RBITC Conjugated

Catalog No: #C92625RBITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VDAC3 Polyclonal Antibody RBITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/283
Conjugates	RBITC
Target Name	VDAC3
Other Names	VDAC3
Accession No.	Uniprot:Q9Y277Gene ID:7419
Uniprot	Q9Y277
GeneID	7419
Calculated MW	31
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only