

DNase I Polyclonal Antibody APC Conjugated

Catalog No: #C92626APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNase I Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/282
Conjugates	APC
Target Name	DNase I
Other Names	DNase I
Accession No.	Uniprot:P24855Gene ID:1773
Uniprot	P24855
GeneID	1773
Calculated MW	29
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only