

## AX2R Polyclonal Antibody FITC Conjugated

Catalog No: #C92630FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	AX2R Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms
Immunogen Type	30-90/193
Conjugates	FITC
Target Name	AX2R
Other Names	AX2R
Accession No.	Uniprot:Q3ZCQ2Gene ID:389289
Uniprot	Q3ZCQ2
GeneID	389289
Calculated MW	22
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=3ug/Test IF=1:50-200

Note: This product is for in vitro research use only