

DR6 Polyclonal Antibody FITC Conjugated

Catalog No: #C92632FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DR6 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/655
Conjugates	FITC
Target Name	DR6
Other Names	DR6
Accession No.	Uniprot:O75509Gene ID:27242
Uniprot	O75509
GeneID	27242
Calculated MW	52/68
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=3ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only