

Corin Polyclonal Antibody Biotin Conjugated

Catalog No: #C92633Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Corin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	165-245/1042
Conjugates	Bio
Target Name	Corin
Other Names	Corin
Accession No.	Uniprot:Q9Y5Q5Gene ID:10699
Uniprot	Q9Y5Q5
GeneID	10699
Calculated MW	115
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only