

P2RX4 Polyclonal Antibody PE Conjugated

Catalog No: #C92635PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	P2RX4 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	238-338/388
Conjugates	PE
Target Name	P2RX4
Other Names	P2RX4
Accession No.	Uniprot:Q99571Gene ID:5025
Uniprot	Q99571
GeneID	5025
Calculated MW	43
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only