

PCNT Polyclonal Antibody FITC Conjugated

Catalog No: #C92638FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | PCNT Polyclonal Antibody FITC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 3101-3200/3336 |
| Conjugates | FITC |
| Target Name | PCNT |
| Other Names | PCNT |
| Accession No. | Uniprot:O95613Gene ID:5116 |
| Uniprot | O95613 |
| GeneID | 5116 |
| Calculated MW | 367 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only