

Factor VII light chain Polyclonal Antibody FITC Conjugated

Catalog No: #C92643FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Factor VII light chain Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	52-150/466
Conjugates	FITC
Target Name	Factor VII light chain
Other Names	Factor VII light chain
Accession No.	Uniprot:P08709Gene ID:2155
Uniprot	P08709
GeneID	2155
Calculated MW	17/28
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only