

## CEND1 Polyclonal Antibody FITC Conjugated

Catalog No: #C92653FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CEND1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	4-100/149
Conjugates	FITC
Target Name	CEND1
Other Names	CEND1
Accession No.	Uniprot:Q8N111Gene ID:51286
Uniprot	Q8N111
GeneID	51286
Calculated MW	15
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only