

DLL3 Polyclonal Antibody Biotin Conjugated

Catalog No: #C92654Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DLL3 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms
Immunogen Type	51-150/618
Conjugates	Bio
Target Name	DLL3
Other Names	DLL3
Accession No.	Uniprot:Q9NYJ7Gene ID:10683
Uniprot	Q9NYJ7
GeneID	10683
Calculated MW	65
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only