

Cyclin B3 Polyclonal Antibody HRP Conjugated

Catalog No: #C92657HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---------------------------------------------------------------------------------|
| Product Name | Cyclin B3 Polyclonal Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | WB IHC |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 1117-1200/1395 |
| Conjugates | HRP |
| Target Name | Cyclin B3 |
| Other Names | Cyclin B3 |
| Accession No. | Uniprot:Q8WWL7Gene ID:85417 |
| Uniprot | Q8WWL7 |
| GeneID | 85417 |
| Calculated MW | 158 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Note: This product is for in vitro research use only