

HHLA2 Polyclonal Antibody PE Conjugated

Catalog No: #C92679PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HHLA2 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu
Immunogen Type	231-330/414
Conjugates	PE
Target Name	HHLA2
Other Names	HHLA2
Accession No.	Uniprot:Q9UM44Gene ID:11148
Uniprot	Q9UM44
GeneID	11148
Calculated MW	45
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only