

BTNL9 Polyclonal Antibody HRP Conjugated

Catalog No: #C92681HRP

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | BTNL9 Polyclonal Antibody HRP Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | WB IHC ICC |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 121-220/535 |
| Conjugates | HRP |
| Target Name | BTNL9 |
| Other Names | BTNL9 |
| Accession No. | Uniprot:Q6UXG8Gene ID:153579 |
| Uniprot | Q6UXG8 |
| GeneID | 153579 |
| Calculated MW | 56 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

Note: This product is for in vitro research use only