

Lin28 Polyclonal Antibody PE Conjugated

Catalog No: #C92683PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lin28 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	75-180/209
Conjugates	PE
Target Name	Lin28
Other Names	Lin28
Accession No.	Uniprot:Q9H9Z2Gene ID:79727
Uniprot	Q9H9Z2
GeneID	79727
Calculated MW	29
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only