

# phospho-ATF2 (Thr53) Polyclonal Antibody HRP Conjugated



Catalog No: #C92684HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	phospho-ATF2 (Thr53) Polyclonal Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	TP(p-T)PT
Conjugates	HRP
Target Name	phospho-ATF2 (Thr53)
Other Names	phospho-ATF2 (Thr53)
Accession No.	Uniprot:P15336Gene ID:1386
Uniprot	P15336
GeneID	1386
Calculated MW	54
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only